

mitogen activated protein kinase 1, MAPK1; MAPK, ELISA Kit

Cat PEK05630

Species

Plant

Full Product Name

Plant mitogen activated protein kinase 1, MAPK1; MAPK ELISA Kit

Product Gene Name

MAPK1 elisa kit

Assay Type

Quantitative Competitive

Detection Range

0.19-12ng/mL

Sensitivity

0.013ng/mL

Samples

Serum, plasma, cell culture supernates, cell lysates, tissue homogenates and other biological fluids.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY