

CASP-like protein 6 (LOC100793096), Recombinant Protein

Cat *RP18466*

Species

Glycine max (Soybean) (*Glycine hispida*)

Full Product Name

Recombinant Glycine max CASP-like protein 6, partial

Product Gene Name

LOC100793096 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

20,187 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

CASP-like protein

NCBI Accession

NP_001239666.1

NCBI GI

358248678

NCBI GenBank Nucleotide

NM_001252737.2

NCBI GeneID

100793096

NCBI Official Full Name

CASP-like protein 1E1

NCBI Official Symbol

LOC100793096

NCBI Protein Information

CASP-like protein 1E1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

CASP-like protein 6 (LOC100793096), Recombinant Protein

Cat *RP18466*

UniProt Gene Name

GmCASPL1E1

UniProt Synonym Gene Names

GmCASPL1E1

UniProt Protein Name

CASP-like protein 1E1

UniProt Primary Accession

C6TBD0

UniProt Related Accession

C6TBD0

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY