

30 kDa cell wall protein, Recombinant Protein

Cat RP00486

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

(Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Arabidopsis thaliana (Mouse-ear cress)
Cell)

Full Product Name

Recombinant Arabidopsis thaliana 30 kDa cell wall protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

NFQRDVEITW GDMRR

Sequence Positions

1-5, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

1,923 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

30 kDa cell wall protein

NCBI Accession

P80853.1

NCBI GI

82592787

NCBI Official Full Name

30 kDa cell wall protein

UniProt Protein Name

30 kDa cell wall protein

UniProt Primary Accession

P80853

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY