

UPF0496 protein At4g34320, Recombinant Protein

Cat *RP00493*

Size *0.02 mg/ 0.1 mg/ 5x0.1 mg*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* UPF0496 protein At4g34320 (At4g34320)

Product Gene Name

AT4G34320 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MGNQTSKKSQ ETSKSVHYT TELRSYAAAC KADTELQSFD TCLQARTSHV ISTLATGVEV RALSFDSLKE
VTQCLLEMNQ EVVKVILDCK KDIWKNQEMF ELVEDYFENS LKTLDFCAAL EKGLRRARDS HLLILVALQQ
FEDESLVQGG NGYKKTLEEL KNFKDAESPF NEDFFKMFQS VYKQQLMLE KLQHRKNKLD KKLKCIHTWR
KLSSIIFVAT FATVLICSVV AAAMAAPPVA AALAAATAVP LGSMGKWIDS LWKNYENALK GQKEVISSMQ
AGTFVAVKDL DNIRVLIERL EIEITGMVKS AEFAVEHNAV KIGIDDIKKK LEVFKKNVEE LGTQADLCSR
DIRRARTVIL QRIIKHPNNA SSST

Sequence Positions

1-374aa; full length protein

Format

Liquid containing glycerol

Host

Cell Free Expression

Molecular Weight

42,054 Da

Storage

Store at -20°C, for extended storage, conserve at -20°C or -80°C. Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Protein Family

UPF0496 protein

NCBI Accession

NP_195158.1

NCBI GI

15235330

NCBI GenBank Nucleotide

NM_119597.3

NCBI GeneID

829582

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

UPF0496 protein At4g34320, Recombinant Protein

Cat *RP00493*

Size *0.02 mg/ 0.1 mg/ 5x0.1 mg*

NCBI Official Full Name

uncharacterized protein

NCBI Official Symbol

AT4G34320

NCBI Official Synonym Symbols

F10M10.90; F10M10_90

NCBI Protein Information

hypothetical protein

UniProt Gene Name

At4g34320

UniProt Protein Name

UPF0496 protein At4g34320

UniProt Entry Name

U496A_ARATH

UniProt Primary Accession

Q9SYZ7

UniProt Related Accession

Q9SYZ7

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
