

CASP-like protein ARALYDRAFT_470341, Recombinant Protein

Cat *RP00499*

Size *0.02 mg/ 0.1 mg/ 5x0.1 mg*

Species

Arabidopsis lyrata subsp. *lyrata* (Lyre-leaved rock-cress)

Full Product Name

Recombinant *Arabidopsis lyrata* subsp. *lyrata* CASP-like protein ARALYDRAFT_470341 (ARALYDRAFT_470341)

Product Gene Name

ARALYDRAFT_470341 recombinant protein

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MVKLTQRLGG LVLRF AAFCA ALGAVIAMIT SRERSSFFVI SLVAKYSDLA AFKYFVIANA IVTVYSFLVL
FLPKESLLWK FVVVLDLMVT MLLTSSLSAA VAVAQVGKRG NANAGWLPIG GQVPRFCDQI TGALIAGLVA
LVLYVFLIF SIHHVDPFL LRKS

Sequence Positions

1-164aa; full length protein

Format

Liquid containing glycerol

Host

Cell Free Expression

Molecular Weight

17,859 Da

Storage

Store at -20°C, for extended storage, conserve at -20°C or -80°C. Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Protein Family

CASP-like protein

NCBI Accession

XP_002889467.1

NCBI GI

297843172

NCBI GenBank Nucleotide

XM_002889421.1

NCBI GeneID

9325531

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

CASP-like protein ARALYDRAFT_470341, Recombinant Protein

Cat *RP00499*

Size *0.02 mg/ 0.1 mg/ 5x0.1 mg*

NCBI Official Full Name

integral membrane family protein

NCBI Official Symbol

ARALYDRAFT_470341

NCBI Protein Information

integral membrane family protein

UniProt Gene Name

ARALYDRAFT_470341

UniProt Synonym Gene Names

AICASPL1C2

UniProt Protein Name

CASP-like protein 1C2

UniProt Entry Name

CSPLE_ARALL

UniProt Primary Accession

D7KCH2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
