

Filament-like plant protein 7 (FPP7), Recombinant Protein

Cat *RP01354*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Filament-like plant protein 7 (FPP7) , partial

Product Gene Name

FPP7 recombinant protein

Product Synonym Gene Name

FPP7

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

101,423 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Filament-like plant protein

NCBI Accession

NP_179917.2

NCBI GI

334184390

NCBI GenBank Nucleotide

NM_127900.4

NCBI GeneID

816868

NCBI Official Full Name

filament-like protein (DUF869)

NCBI Official Symbol

AT2G23360

NCBI Official Synonym Symbols

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Filament-like plant protein 7 (FPP7), Recombinant Protein



Cat *RP01354*

F26B6.1; F26B6_1

NCBI Protein Information

filament-like protein (DUF869)

UniProt Gene Name

FPP7

UniProt Synonym Gene Names

AtFPP7

UniProt Protein Name

Filament-like plant protein 7

UniProt Primary Accession

Q9SLN1

UniProt Related Accession

Q9SLN1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY