

Bauhinia purpurea Lectin (BPL/BPA) - MagneZoom™

Cat PL00421

Size 1 mL/5 mL

Sugar Specificity

Gal β 3GalNAc

Description

Affinity purified Bauhinia lectin (BPL/BPA) is fixed on aminated paramagnetic beads with EDC and stabilized in buffer to maintain maximum activity. These beads are ideal for iP applications to pull down glycoproteins.

Abbreviation

BPL/BPA

Material Source

Camel's foot tree seeds

Conjugate

MagneZoom™

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

Gal β 3GalNAc

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY