

# Dolichos biflorus Lectin (DBA) - MagneZoom™

Cat PL00434

Size 1 mL/5 mL

---

## Sugar Specificity

$\alpha$ GalNAc

## Description

Affinity purified biichos lectins (DBA) are fixed on aminated paramagnetic beads with EDC and stabilized in the buffer to maintain maximum activity. These beads are ideal for iP applications to pull down glycoproteins.

## Abbreviation

DBA

## Material Source

Horse gram seeds

## Conjugate

MagneZoom™

## Shelf Life

1 year

## Blood Group Specificity

Non-specific

## Sugar Specificity

$\alpha$ GalNAc

## Mitogenic Activity

No

## Format

Liquid

## Storage

2-8°C

## Hazardous Shipping

Non-hazardous

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**