

# Glycine max Lectin (SBA) - Coated Microplates

Cat *PL00445*

Size *5 Plate(s)/20 Plate(s)*

---

## Sugar Specificity

### Description

The affinity purified glycine max(soybean) lectin (SBA) is fixed on a microplate for enzyme-linked immunosorbent assay. Unoccupied charged groups are blocked by BSA to reduce non-specific interactions and reduce background.

### Shelf Life

2 years unopened, 2 weeks opened

### Storage

2-8°C

### Hazardous Shipping

Non-hazardous

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---