

Glycine max Lectin (SBA) - MagneZoom™

Cat PL00453

Size 1 mL/5 mL

Sugar Specificity

α - β GalNAc

Description

Affinity purified glycine max(soybean) lectin (SBA) is fixed on aminated paramagnetic beads using EDC and stabilized in buffer to maintain maximum activity. These beads are ideal for iP applications to pull down glycoproteins.

Abbreviation

SBA

Material Source

Soybean seeds

Conjugate

MagneZoom™

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

α - β GalNAc

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY