

Griffonia simplicifolia Lectin (GSL I) - MagneZoom™

Cat PL00470

Size 2 mL/5 mL

Sugar Specificity

α Gal, α GalNAc

Description

Affinity-purified Griffonia simplicifolia lectins (GSLI) are immobilized on aminated paramagnetic beads with EDC and stabilized in buffers to maintain maximum activity. These beads are ideal for iP applications to pull down glycoproteins.

Abbreviation

GSL I

Material Source

Griffonia simplicifolia seeds

Conjugate

MagneZoom™

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

α Gal, α GalNAc

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY