Sambucus nigra Lectin (SNA/EBL II) - Separopore® 4B



Cat PL00571

Size 10 mL/25 mL

Sugar Specificity

GalNAc, Galactose

Description

After CNBr coupling, affinity purified black ginseng lectin (SNA/EBL II) is fixed on Separopore®4B. The final product is stabilized in the buffer to maintain maximum activity. The gel can be filled with 1 ml or 5 ml.

Abbreviation

SNA/EBL II

Material Source

Elderberry bark

Conjugate

Separopore® 4B

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

GalNAc, Galactose

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Address: SUITE 209, 17 Ramsey Road, Shirley, NY 11967

E-mail: info@cd-biosci.com

Tel: 1-631-637-0420

https://www.cd-biosciences.com/plant-protein/