

Phaseolus vulgaris Lectin (PHA-E) - Separopore® 4B

Cat PL00587

Size 2 mL/5 mL/10 mL

Sugar Specificity

Gal β 4GlcNAc α 2Man α 6 (GlcNAc β 4) (GlcNAc β 4Man α 3) Man β 4

Description

Affinity purified bean lectin (PHA-E) is fixed on Separopore®4B after CNBr coupling. The final product is stabilized in the buffer to maintain maximum activity. The gel can be filled with 1 ml or 5 ml.

Abbreviation

PHA-E

Material Source

Red kidney bean seeds

Conjugate

Separopore® 4B

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

Gal β 4GlcNAc α 2Man α 6 (GlcNAc β 4) (GlcNAc β 4Man α 3) Man β 4

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY