

Phaseolus vulgaris Lectin (PHA-E) - Separopore® 6B SuperTrap™ Cartridge



Cat PL00590

Size 3 mL/6 mL/5 mL/10 mL

Sugar Specificity

Description

The crosslinked 6B Separopore® is properly coupled to the Phaseolus vulgaris lectin (phae), providing a fast-flowing cartridge that can capture the affinity-bound material and retain it during washing, then elute the affinity-purified material.

Purity

High Purity Grade

Shelf Life

1 year

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY