

Phaseolus vulgaris Lectin (PHA-E+L) - Separopore® 6B SuperTrap™ Cartridge

Cat PL00607

Size 3 mL/6 mL/5 mL/10 mL

Sugar Specificity

Description

The crosslinked 6B Separopore® is properly coupled to the Phaseolus vulgaris lectin (PHA-E+L), providing a fast-flowing cartridge that can capture the affinium-bound material and retain it during washing, then elute the affinium-purified material.

Purity

High Purity Grade

Shelf Life

1 year

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY