

# Terpenoid synthase 15 (TPS15), Recombinant Protein

Cat      *RP05103*

---

## Species

*Arabidopsis thaliana* (Mouse-ear cress)

## Full Product Name

Recombinant *Arabidopsis thaliana* Terpenoid synthase 15 (TPS15) , partial

## Product Gene Name

TPS15 recombinant protein

## Product Synonym Gene Name

TPS15

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

69,668 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

Terpenoid synthase

## NCBI Accession #

NP\_189564.2

## NCBI GI #

186510561

## NCBI GenBank Nucleotide #

NM\_113843.3

## NCBI GeneID

822572

## NCBI Official Full Name

Terpenoid cyclases/Protein prenyltransferases superfamily protein

## NCBI Official Symbol

AT3G29190

## NCBI Protein Information

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# Terpenoid synthase 15 (TPS15), Recombinant Protein

Cat      *RP05103*

---

Terpenoid cyclases/Protein prenyltransferases superfamily protein

## UniProt Gene Name

TPS15

## UniProt Synonym Gene Names

AtTPS15

## UniProt Protein Name

Terpenoid synthase 15

## UniProt Primary Accession #

Q9LS76

## UniProt Related Accession #

Q9LS76

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**