

# B3 domain-containing protein REM10 (REM10), Recombinant Protein

Cat *RP05141*

---

## Species

*Arabidopsis thaliana* (Mouse-ear cress)

## Full Product Name

Recombinant *Arabidopsis thaliana* B3 domain-containing protein REM10 (REM10) , partial

## Product Gene Name

REM10 recombinant protein

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

63,258 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

B3 domain-containing protein

## NCBI Accession #

NP\_180046.2

## NCBI GI #

42569291

## NCBI GenBank Nucleotide #

NM\_128031.4

## NCBI GenelD

817006

## NCBI Official Full Name

Transcriptional factor B3 family protein

## NCBI Official Symbol

AT2G24700

## NCBI Official Synonym Symbols

F27A10.1; F27A10\_1

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# B3 domain-containing protein REM10 (REM10), Recombinant Protein

Cat *RP05141*

---

## NCBI Protein Information

Transcriptional factor B3 family protein

## UniProt Gene Name

REM10

## UniProt Protein Name

B3 domain-containing protein REM10

## UniProt Synonym Protein Names

Protein REPRODUCTIVE MERISTEM 10

## UniProt Primary Accession #

P0DH85

## UniProt Secondary Accession #

Q9SHV3

## UniProt Related Accession #

P0DH85

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**