

# Defensin-like protein 71 (LCR84), Recombinant Protein

Cat      *RP05158*

Size      0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg  
(Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

Species (Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg  
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-  
Arabidopsis thaliana (Mouse-ear cress)  
Cell)

## Full Product Name

Recombinant Arabidopsis thaliana Defensin-like protein 71 (LCR84)

## Product Gene Name

LCR84 recombinant protein

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Sequence

LPSSVINGFG YDYCIVECSI TFLDDACKPL CIDRGYSDGG CIGLSPHFKC CCKK

## Sequence Positions

23-76, Full length protein

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

8,227 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week.  
Repeated freezing and thawing is not recommended.

## Protein Family

Defensin-like protein

## NCBI Accession #

NP\_001031355.1

## NCBI GI #

79322334

## NCBI GenBank Nucleotide #

NM\_001036278.3

## NCBI GeneID

3768590

## NCBI Official Full Name

low-molecular-weight cysteine-rich 84

## NCBI Official Symbol

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

# Defensin-like protein 71 (LCR84), Recombinant Protein

Cat      *RP05158*

Size      0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg  
(Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

LCR84      (Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg  
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-  
Cell)

## NCBI Official Synonym Symbols

low-molecular-weight cysteine-rich 84

## NCBI Protein Information

low-molecular-weight cysteine-rich 84

## NCBI Summary

Encodes a member of a family of small, secreted, cysteine rich protein with sequence similarity to the PCP (pollen coat protein) gene family.

## UniProt Gene Name

LCR84

## UniProt Synonym Gene Names

Protein LCR84

## UniProt Protein Name

Defensin-like protein 71

## UniProt Synonym Protein Names

Low-molecular-weight cysteine-rich protein 84; Protein LCR84

## UniProt Primary Accession #

P82793

## UniProt Related Accession #

P82793

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**