

Autophagy-related protein 8b (ATG8B), Recombinant Protein

Cat RP05167

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

Species (Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana Autophagy-related protein 8b (ATG8B)

Product Gene Name

ATG8B recombinant protein

Product Synonym Gene Name

ATG8B

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MEKNSFKLSN PLEMRMAEST RIRAKYPERV PVIVEKAGQS DVPDIDKKKY LVPADLTIGQ FVYVVRKRIK
LGAEKAIFVF VKNTLPPTAA LMSAIYEEHK DEDGFLYMTY SGENTFG

Sequence Positions

1-117, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

13,880 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Autophagy-related protein

NCBI Accession

NP_001329707.1

NCBI GI

1063720985

NCBI GenBank Nucleotide

NM_001340494.1

NCBI GeneID

825794

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Autophagy-related protein 8b (ATG8B), Recombinant Protein

Cat *RP05167*

Size 0.02 mg (E-Coli)/ 0.1 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 1 mg (E-Coli)/ 0.02 mg

(Mammalian-Cell)/ 0.1 mg (Baculovirus)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Cell)
Ubiquitin-like superfamily protein

NCBI Official Symbol

ATG8B

NCBI Official Synonym Symbols

autophagy 8b; F4H6.14; F4H6_14

NCBI Protein Information

Ubiquitin-like superfamily protein

UniProt Gene Name

ATG8B

UniProt Synonym Gene Names

APG8B; AtAPG8b; Protein autophagy 8b

UniProt Protein Name

Autophagy-related protein 8b

UniProt Synonym Protein Names

Autophagy-related ubiquitin-like modifier ATG8b; AtAPG8b; Protein autophagy 8b

UniProt Primary Accession

Q9XEB5

UniProt Related Accession

Q9XEB5

UniProt Comments

Ubiquitin-like modifier involved in autophagosomes formation. May mediate the delivery of the autophagosomes to the vacuole via the microtubule cytoskeleton.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY