

# F-box/LRR-repeat protein 10 (FBL10), Recombinant Protein

Cat RP05306

---

## Species

Arabidopsis thaliana (Mouse-ear cress)

## Full Product Name

Recombinant Arabidopsis thaliana F-box/LRR-repeat protein 10 (FBL10) , partial

## Product Gene Name

FBL10 recombinant protein

## Product Synonym Gene Name

FBL10

## Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

## Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

## Host

E Coli or Yeast or Baculovirus or Mammalian Cell

## Molecular Weight

72,282 Da

## Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

## Protein Family

F-box/LRR-repeat protein

## NCBI Accession #

NP\_565400.1

## NCBI GI #

18398283

## NCBI GenBank Nucleotide #

NM\_127255.4

## NCBI GeneID

816205

## NCBI Official Full Name

F-box/RNI-like superfamily protein

## NCBI Official Symbol

AT2G17020

## NCBI Official Synonym Symbols

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# F-box/LRR-repeat protein 10 (FBL10), Recombinant Protein

Cat *RP05306*

---

F6P23.18; F6P23\_18

## NCBI Protein Information

F-box/RNI-like superfamily protein

## UniProt Gene Name

FBL10

## UniProt Protein Name

F-box/LRR-repeat protein 10

## UniProt Primary Accession #

Q9SDA8

## UniProt Secondary Accession #

Q8LF00

## UniProt Related Accession #

Q9SDA8

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**