

Tubby-like F-box protein 6 (TULP6), Recombinant Protein

Cat *RP05377*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana Tubby-like F-box protein 6 (TULP6)

Product Gene Name

TULP6 recombinant protein

Product Synonym Gene Name

TULP6

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MDLSLLLPRK SWEISGGADK ISRVAMSLKN IVKNKYKAIG RRGSRSHIAPE GSSVSSSLST NEGLNQSIWV
DLPPPELLLDI IQRIESEQSL WPGRRDVVAC ASVCKSWREM TKEVVVKVPEL SGLITFPISL RQPGPRDAPI
QCFIKRERAT GIYRLYLGLS PALSGDKSKL LLSAKRVRRR TGAEFVVSLS GNDFSRSSSN YIGKLRSNFL
GTKFTVYENQ PPPFNRLKPP SMQVSPWVSS SSSSYNIASI LYELNVLRTG GPRRMQCIMH SIPISAIQEG
GKIQSPTFT NQGKKKKKPL MDFCSGNLGG ESVIKEPLIL KNKSPRWHEQ LQCWCLNFKG RVTVASVKNF
QLVAAAAEAG KMNPIPEEEQ DRVILQFGKI GKDIFTMDYR YPISAFQAFI ICLSSFDTKP VCE

Sequence Positions

1-413, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

46,065 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Tubby-like F-box protein

NCBI Accession

NP_001322089.1

NCBI GI

1063693900

NCBI GenBank Nucleotide

NM_001333268.1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Tubby-like F-box protein 6 (TULP6), Recombinant Protein

Cat *RP05377*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg
(Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg
(Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-
Cell)

NCBI GeneID

841130

NCBI Official Full Name

tubby like protein 6

NCBI Official Symbol

TLP6

NCBI Official Synonym Symbols

AtTLP6; F8G22.1; F8G22_1; tubby like protein 6

NCBI Protein Information

tubby like protein 6

NCBI Summary

Member of TLP family

UniProt Gene Name

TULP6

UniProt Synonym Gene Names

AtTLP6

UniProt Protein Name

Tubby-like F-box protein 6

UniProt Primary Accession

Q0WPY0

UniProt Secondary Accession

Q9C6B4

UniProt Related Accession

Q0WPY0

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY