

Exocyst complex component 5 (SEC10), Recombinant Protein

Cat *RP05420*

Species

Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant *Arabidopsis thaliana* Exocyst complex component 5 (SEC10) , partial

Product Gene Name

SEC10 recombinant protein

Product Synonym Gene Name

SEC10

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Exocyst complex component

NCBI Accession

Q8RVQ5.2

NCBI GI

152013638

NCBI GenelD

831112

NCBI Official Symbol

SEC10

NCBI Official Synonym Symbols

exocyst complex component sec10

NCBI Protein Information

exocyst complex component sec10

UniProt Primary Accession

Q8RVQ5

UniProt Related Accession

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Exocyst complex component 5 (SEC10), Recombinant Protein

Cat *RP05420*

X5JA13

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY