

Putative F-box only protein 11 (FBX11), Recombinant Protein

Cat *RP05494*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)
Species Arabidopsis thaliana (Mouse-ear cress)

Full Product Name

Recombinant Arabidopsis thaliana Putative F-box only protein 11 (FBX11)

Product Gene Name

FBX11 recombinant protein

Product Synonym Gene Name

FBX11

Purity

Greater or equal to 85% purity as determined by SDS-PAGE. (lot specific)

Sequence

MVSVNLPWEL VEEILCRVPP QSLVKFRTVC KQWNSLFDDN KFVNDHFVQS QPQLIFRTES KIYSVAVNFK
GPRIEVHELP LAIPGLKSEM PIRLHDYVDC DGLLFCTSYF NGVLIWNPWL RQTRFFPSIH RYPVTYDIGY
DNKKQYKMLD YYKCEGDSSI KSIIYEIGLY SKKVKDLESD STWSIFQSNS VSLNGTLYWA GVDVNNGIFI
RSFSFSTEKP RTFCSLPFMY DDDNVLALAV FRKDRLSLN LCNKTSEIKI WLTKNNINDR EVGLEEDVWW
INLMTVLIPN FPKFSFYWYN RPDLYFLDN DDAKRLVICC YDETDQAYIY IVKGDMFKKI KIDQYEVLS
YSRPHFHTYI PSLVRPSRVK EDNKN

Sequence Positions

1-375, Full length protein

Format

Lyophilized or liquid (Format to be determined during the manufacturing process)

Host

E Coli or Yeast or Baculovirus or Mammalian Cell

Molecular Weight

44,208 Da

Storage

Store at -20°C. For long-term storage, store at -20°C or -80°C. Store working aliquots at 4°C for up to one week. Repeated freezing and thawing is not recommended.

Protein Family

Putative F-box only protein

NCBI Accession

NP_564587.1

NCBI GI

18403471

NCBI GenBank Nucleotide

NM_104010.1

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Putative F-box only protein 11 (FBX11), Recombinant Protein

Cat *RP05494*

Size 0.02 mg (E-Coli)/ 0.02 mg (Yeast)/ 0.1 mg (E-Coli)/ 0.1 mg (Yeast)/ 0.02 mg (Baculovirus)/ 0.02 mg (Mammalian-Cell)/ 0.1

mg (Baculovirus)/ 1 mg (E-Coli)/ 1 mg (Yeast)/ 0.1 mg (Mammalian-Cell)/ 1 mg (Baculovirus)/ 0.5 mg (Mammalian-Cell)

NCBI GeneID

841555

NCBI Official Full Name

F-box and associated interaction domains-containing protein

NCBI Official Symbol

AT1G51320

NCBI Official Synonym Symbols

F11M15.18; F11M15_18

NCBI Protein Information

F-box and associated interaction domains-containing protein

UniProt Gene Name

FBX11

UniProt Protein Name

Putative F-box only protein 11

UniProt Primary Accession

Q9SYD4

UniProt Related Accession

Q9SYD4

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY