

DYKDDDDK (mouse monoclonal) (clone F-tag-01)

Cat PA00003

Size 0.1 mg

Host

Mouse

Clonality

Monoclonal

Confirmed reactivity

DYKDDDDK (Sigma FLAG) epitope tag

Immunogen

KLH-conjugated synthetic peptide: DYKDDDDK (Sigma FLAG)

Host

Mouse

Clonality

Monoclonal

Purity

Immunoglobulin purified on Protein A in PBS pH 7.4. Contains 15 mM sodium azide.

Format

Liquid at 1 mg/ml.

Storage

Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application

Immunocyto chemistry (ICC), Western blot (WB)

Recommended dilution

1 µg/ml (WB)

Confirmed reactivity

DYKDDDDK (Sigma FLAG) epitope tag

Additional information

This antibody recognizes DYKDDDDK (Sigma FLAG) tagged proteins

Description

Mouse anti-DyKDDDDK is a primary antibody (Sigma FLAG) that binds to DYKDDDDK. The label's small size and high hydrophilicity reduce the likelihood of interference with its expression, proteolytic maturation, antigenicity, localization, and function.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY