

DYKDDDDK (polyclonal)

Cat PA00004

Size 50 µg

Host

Rabbit

Clonality

Polyclonal

Confirmed reactivity

DYKDDDDK epitope tag (Sigma FLAG), fused to proteins in plant cells

Immunogen

KLH-conjugated synthetic peptide: DYKDDDDK (Sigma FLAG)

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified serum, in PBS pH 7.4.

Format

Lyophilized

Reconstitution

For reconstitution add 50 µl, of sterile water.

Storage

Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application

Western blot (WB)

Recommended dilution

1 : 1000 - 1: 5000 (WB)

Expected | apparent MW

Depends upon fusion protein partner

Confirmed reactivity

DYKDDDDK epitope tag (Sigma FLAG), fused to proteins in plant cells

Not reactive in

No confirmed exceptions from predicted reactivity are currently known

Description

Rabbit anti-DYKDDDDK is a primary antibody combined with DYKDDDDK (Sigma FLAG).

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY