

His-tag | 3xHis (polyclonal) - Biotin conjugated 40 µg

Cat PA00083

Size 40 µg

Host

Rabbit

Clonality

Polyclonal

Confirmed reactivity

3xHis tagged proteins

Immunogen

KLH-conjugated synthetic peptide 3xHis

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified serum in PBS pH 7.4, conjugated to biotin.

Format

Liquid

Storage

Store at 4°C for 12-18 months. A preservative may be added for long time storage up to 2 years.

Application

Western blot (WB)

Recommended dilution

1: 1000 - 1 : 5000 (WB)

Expected | apparent MW

3xHis tagged proteins

Confirmed reactivity

3xHis tagged proteins

Description

His-Tag is a multi-histidine tag that is introduced by 3 or 6 histidine residues to the N or c terminal of a protein. The multi-histidine label can be used for recombinant protein detection using specific antibodies that do not bind to any dye or enzyme.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY