

# His-tag | 6xHis, conjugated to Europium (Clone AD1,1,10)

Cat PA00098

Size 50 µg (1 mg/ml)

---

## Host

Mouse

## Clonality

Monoclonal

## Confirmed reactivity

## Immunogen

Polyhistidine (HHHHHH) tagged polypeptide

## Host

Mouse

## Clonality

Monoclonal

## Purity

Immunogen affinity purified in 50 mM tris (pH 7.8), 0.9% NaCl, 0.05% sodium azide.

## Format

Liquid

## Storage

Store at 2-8°C. The shelf life is 1 year from date of receipt. Prepare working dilution prior to use and then discard.

## Application

High Throughput Screening (HTS)

## Recommended dilution

The optimal working dilution should be determined by the investigator

## Additional information

Antibody is provided in: 50 mM Tris (pH 7,8), 0,9% Sodium Chloride, 0,05% Sodium Azide

## Description

Mouse anti-6X histidine is a primary antibody that binds the 6x histidine label and is directly coupled to europium. This is due to the advantage of shortening the assay time by not requiring the use of 6x histidine label secondary antibodies.

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY