

BC2-tag - DyLight 488 conjugated (40 µg)

Cat PA00133

Size 40 µg

Host

Rabbit

Clonality

Polyclonal

Confirmed reactivity

Protein of interest overexpressed with BC2-tag

Immunogen

KLH-conjugated BC2 tag PDRKAAVSHWQQ derived from beta catenin.

Host

Rabbit

Clonality

Polyclonal

Purity

Antigen affinity purified serum

Format

Liquid in PBS pH 7,4.

Storage

Store at 4°C for 12-18 months, A preservative may be added for long time storage up to 2 years. Spin briefly the tube before use.

Application

Immunofluorescence (IF), Western blot (WB)

Recommended dilution

To be determined by end user.

Confirmed reactivity

Protein of interest overexpressed with BC2-tag

Additional information

This antibody is suitable for detection of recombinant proteins with BC2 tag,

Description

The BC2 tag is a short, inert peptide tag used in protein studies. DyLight 488 Amax = 493 nm, Emax = 519 nm. DyLight is a registered trademark of Thermofisher Inc. and its subsidiaries.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY