

# BC2-tag - DyLight 650 conjugated (40 µg)

Cat PA00135

Size 40 µg

---

## Host

Rabbit

## Clonality

Polyclonal

## Confirmed reactivity

Protein of interest overexpressed with BC2-tag

## Immunogen

KLH-conjugated BC2 tag PDRKAAVSHWQQ derived from beta catenin.

## Host

Rabbit

## Clonality

Polyclonal

## Purity

Antigen affinity purified serum

## Format

Liquid in PBS pH 7,4.

## Storage

Store at 4°C for 12-18 months, A preservative may be added for long time storage up to 2 years. Spin briefly the tube before use.

## Application

Immunofluorescence (IF), Western blot (WB)

## Recommended dilution

To be determined by end user.

## Confirmed reactivity

Protein of interest overexpressed with BC2-tag

## Additional information

This antibody is suitable for detection of recombinant proteins with BC2 tag,

## Description

The BC2 tag is a short, inert peptide tag used in protein studies. DyLight 650 Amax = 652 nm, Emax = 672 nm. IgGDyLight is a registered trademark of Thermofisher Inc. and its subsidiaries.

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**