

Lotus tetragonolobus Lectin (LTL) - Separopore® 4B

Cat PL00030

Size 2 mL/5 mL

Sugar Specificity

α Fucose

Description

After CNBr coupling, affinity-purified Lotus Tetragonal Lectin (LTL) is immobilized on Separopore® 4B. The final product is stabilized in buffer to maintain maximum activity. The gel can be filled with 1 mL or 5 mL.

Abbreviation

LTL

Material Source

Winged pea

Conjugate

Separopore® 4B

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

α Fucose

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY