

# GST conjugated to Europium 1024

Cat PA00166

Size 50 µg (1 mg/ml)



---

## Host

Goat

## Clonality

Polyclonal

## Confirmed reactivity

## Immunogen

Glutathione-S-transferase from *S. japonicum*

## Host

Goat

## Clonality

Polyclonal

## Purity

Immunogen affinity purified using solid phase GST in 50 mM tris (pH 7.8), 0.9% NaCl, 0.05% sodium azide.

## Format

Liquid

## Storage

Store at 2-8°C. Shelf life: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

## Application

High Throughput Screening (HTS)

## Recommended dilution

The optimal working dilution should be determined by the investigator

## Additional information

Antibody is provided in: 50 mM Tris (pH 7,8), 0,9% Sodium Chloride, 0,05% Sodium Azide

## Description

Goat anti-GST is a primary antibody that binds to GST and is directly coupled to Europium 1024.

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**